

Product datasheet

Anti-CYP20A1 antibody ab113006

1 图像

概述

| | |
|-------|---|
| 产品名称 | Anti-CYP20A1抗体 |
| 描述 | 兔多克隆抗体to CYP20A1 |
| 经测试应用 | 适用于: WB |
| 种属反应性 | 与反应: Human |
| 免疫原 | Synthetic peptide conjugated to KLH, corresponding to a region within internal sequence amino acids 228-258 of Human CYP20A1 (NP_803882.1). |
| 阳性对照 | 293 cell lysate nontransfected and transiently transfected with the CYP20A1 gene. |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term. |
| 存储溶液 | Preservative: 0.09% Sodium azide Constituent: 99% PBS |
| 纯度 | Ammonium Sulphate Precipitation |
| 纯化说明 | ab113006 was prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS. |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab113006** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----|------|---|
| WB | | 1/50 - 1/100. Predicted molecular weight: 52 kDa. |

靶标


序列相似性

Belongs to the cytochrome P450 family.

细胞定位

Membrane.

图片



All lanes : Anti-CYP20A1 antibody
(ab113006) at 1/50 dilution

Lane 1 : non-transfected 293 cell lysate

Lane 2 : CYP20A1 transiently transfected 293
cell lysate

Lysates/proteins at 2 μ g per lane.

Predicted band size : 52 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors