

Product datasheet

Anti-CYP20A1 antibody ab113006

1 图像

概述

产品名称	Anti-CYP20A1抗体
描述	兔多克隆抗体to CYP20A1
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region within internal sequence amino acids 228-258 of Human CYP20A1 (NP_803882.1).
阳性对照	293 cell lysate nontransfected and transiently transfected with the CYP20A1 gene.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium azide Constituent: 99% PBS
纯度	Ammonium Sulphate Precipitation
纯化说明	ab113006 was prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab113006** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/50 - 1/100. Predicted molecular weight: 52 kDa.

靶标

序列相似性	Belongs to the cytochrome P450 family.
细胞定位	Membrane.

图片



All lanes : Anti-CYP20A1 antibody
(ab113006) at 1/50 dilution

Lane 1 : non-transfected 293 cell lysate

Lane 2 : CYP20A1 transiently transfected 293
cell lysate

Lysates/proteins at 2 μ g per lane.

Predicted band size: 52 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors