


Product datasheet

Anti-A RAF antibody ab113271

1 图像

概述

产品名称	Anti-A RAF抗体
描述	兔多克隆抗体to A RAF
宿主	Rabbit
经测试应用	适用于: IP 不适用于: WB
种属反应性	与反应: Human 预测可用于: Sheep, Rabbit, Cow, Dog, Pig, Chimpanzee, Rhesus monkey, Gorilla, Orangutan, Bat 
免疫原	Synthetic peptide, corresponding to a region within N terminal amino acids 1-50 of Human A RAF (NP_001645.1).
阳性对照	HeLa whole cell lysates.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: 99% Tris citrate/phosphate, 0.09% Sodium chloride Note: pH 7 to 8.
纯度	Immunogen affinity purified
纯化说明	ab113271 was affinity purified using an epitope specific to A RAF immobilized on solid support.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab113271** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
IP		Use at 2-10 µg/mg of lysate.

应用说明 Is unsuitable for WB.

靶标

功能	Involved in the transduction of mitogenic signals from the cell membrane to the nucleus.
组织特异性	Predominantly in urogenital tissues.
序列相似性	Belongs to the protein kinase superfamily. TKL Ser/Thr protein kinase family. RAF subfamily. Contains 1 phorbol-ester/DAG-type zinc finger. Contains 1 protein kinase domain. Contains 1 RBD (Ras-binding) domain.

图片



Detection of A RAF by Western Blot of Immunoprecipitate.
anti-A RAF at 1µg/ml staining A RAF in HeLa whole cell lysate immunoprecipitated using ab113271 at 6µg/mg lysate (1 mg/IP; 20% of IP loaded/lane).
Detection: Chemiluminescence with exposure time of 3 seconds.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors