


Product datasheet

Anti-A RAF antibody ab113271

1 图像

概述

产品名称	Anti-A RAF抗体
描述	兔多克隆抗体to A RAF
经测试应用	适用于: IP 不适用于: WB
种属反应性	与反应: Human 预测可用于: Sheep, Rabbit, Cow, Dog, Pig, Chimpanzee, Rhesus monkey, Gorilla, Orangutan, Bat 
免疫原	Synthetic peptide, corresponding to a region within N terminal amino acids 1-50 of Human A RAF (NP_001645.1).
阳性对照	HeLa whole cell lysates.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: 99% Tris citrate/phosphate, 0.09% Sodium chloride Note: pH 7 to 8.
纯度	Immunogen affinity purified
纯化说明	ab113271 was affinity purified using an epitope specific to A RAF immobilized on solid support.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab113271** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IP		Use at 2-10 µg/mg of lysate.

靶标

功能

Involved in the transduction of mitogenic signals from the cell membrane to the nucleus.

组织特异性

Predominantly in urogenital tissues.

序列相似性


Belongs to the protein kinase superfamily. TKL Ser/Thr protein kinase family. RAF subfamily.

Contains 1 phorbol-ester/DAG-type zinc finger.

Contains 1 protein kinase domain.

Contains 1 RBD (Ras-binding) domain.

图片



Detection of A RAF by Western Blot of Immunoprecipitate.

anti-A RAF at 1 μ g/ml staining A RAF in HeLa whole cell lysate immunoprecipitated using ab113271 at 6 μ g/mg lysate (1 mg/IP; 20% of IP loaded/lane).

Detection: Chemiluminescence with exposure time of 3 seconds.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors