

Product datasheet

Anti-PKIB antibody ab113616

★★★★☆ 1 Abreviews 1 图像

概述

产品名称	Anti-PKIB抗体
描述	兔多克隆抗体to PKIB
经测试应用	适用于: WB
种属反应性	与反应: Rat 预测可用于: Mouse, Rabbit, Goat, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Zebrafish
免疫原	Synthetic peptide corresponding to a region within internal sequence amino acids 40-89 (TDVESVISSF ASSARAGR RN ALPDIQSSLA TGGSPDLALK LEALAVKEDA) of Rat PKIB (NP_001070021). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Rat Kidney Lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: 98% PBS, 2% Sucrose
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab113616** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB	★★★★☆	Use a concentration of 1 µg/ml. Predicted molecular weight: 8 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.
----	-------	---

靶标


功能

Extremely potent competitive inhibitor of cAMP-dependent protein kinase activity, this protein interacts with the catalytic subunit of the enzyme after the cAMP-induced dissociation of its regulatory chains.

序列相似性

Belongs to the PKI family.

图片



Anti-PKIB antibody (ab113616) at 1 μ g/ml +
Rat kidney lysate at 10 μ g

Predicted band size : 8 kDa

10-20% SDS PAGE

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors