

Product datasheet

Anti-Myeloid Precursor Cells antibody [OX-82] (Phycoerythrin)  
ab114040

1 References 1 图像

概述

产品名称	Anti-Myeloid Precursor Cells抗体[OX-82] (Phycoerythrin)
描述	小鼠单克隆抗体[OX-82] to Myeloid Precursor Cells (Phycoerythrin)
偶联物	Phycoerythrin. Ex: 488nm, Em: 575nm
特异性	ab114040 recognizes a 35 kDa antigen found on myeloid cells and stromal elements from a variety of tissues in the adult rat.
经测试应用	适用于: WB, IHC-Fr, Flow Cyt
种属反应性	与反应: Rat
免疫原	Anemic rat bone marrow
阳性对照	Wistar Rat bone marrow cells

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.02% Sodium azide Constituents: 0.5% BSA, 99.48% PBS
纯度	Protein G purified
克隆	单克隆
克隆编号	OX-82
同种型	IgG1

应用


Our [Abpromise guarantee](#) covers the use of **ab114040** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use at an assay dependent dilution. Predicted molecular weight: 35 kDa.

应用	Ab评论	说明
IHC-Fr		Use at an assay dependent dilution.
Flow Cyt		Use 0.5µg for 10 <sup>6</sup> cells.

## 图片



Flow cytometry image staining Wistar Rat bone marrow using ab114040 at 0.5µg per 10<sup>6</sup> cells. Percentage of cells stained above control: 21.5%.

Flow Cytometry - Anti-Myeloid Precursor Cells antibody [OX-82] (Phycoerythrin) (ab114040)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors