

Product datasheet

# Anti-CNNM4 antibody ab116131

1 图像

概述

产品名称	Anti-CNNM4抗体
描述	兔多克隆抗体to CNNM4
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region within internal sequence amino acids 500-530 of Human CNNM4 (NP_064569.3).
阳性对照	CEM cell line lysates.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium azide Constituent: 99% PBS
纯度	Immunogen affinity purified
纯化说明	ab116131 is purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab116131** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/100 - 1/500. Predicted molecular weight: 87 kDa.

靶标

功能 Probable metal transporter. The interaction with the metal ion chaperone COX11 suggests that it

may play a role in sensory neuron functions (By similarity). May play a role in biomineralization and retinal function.

**组织特异性**

Widely expressed. Highly expressed in heart.

**疾病相关**

Defects in CNNM4 are the cause of Jalili syndrome (JALIS) [MIM:217080]. A syndrome characterized by the association of cone-rod dystrophy and amelogenesis imperfecta.


**序列相似性**

Belongs to the ACDP family.  
Contains 2 CBS domains.  
Contains 1 DUF21 domain.

**细胞定位**

Cell membrane.

**Anti-CNNM4 antibody 图像**



Anti-CNNM4 antibody (ab116131) at 1/100 dilution + CEM cell line lysate at 35 µg

**Predicted band size : 87 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors