

Product datasheet

Anti-Zfp472 antibody ab116228

1 图像

概述

| | |
|-------|---|
| 产品名称 | Anti-Zfp472抗体 |
| 描述 | 兔多克隆抗体to Zfp472 |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: WB |
| 种属反应性 | 与反应: Mouse 预测可用于: Rat, Horse, Cow, Cat, Dog, Pig |
| 免疫原 | Synthetic peptide corresponding to a region within N terminal amino acids 20-69 (TLGEWTMLDS SQKKLYRDVM KETFLNLISI GKTEEDTEE EYQNPKGNLR) of Mouse Zfp472 (NP_694703.2). Run BLAST with ExPASy Run BLAST with NCBI |
| 阳性对照 | Mouse brain lysate. |

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Constituents: 97% PBS, 2% Sucrose |
| 纯度 | Immunogen affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用


Our [Abpromise guarantee](#) covers the use of **ab116228** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----|------|--|
| WB | | Use a concentration of 1 µg/ml. Predicted molecular weight: 62 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T. |

靶标

图片



Anti-Zfp472 antibody (ab116228) at 1 μ g/ml +
Mouse brain lysate at 10 μ g

Predicted band size: 62 kDa

Gel concentration 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors