

Product datasheet

Anti-Ephrin A2 antibody [2C6] ab117941

1 图像

概述

产品名称	Anti-Ephrin A2抗体[2C6]
描述	小鼠单克隆抗体[2C6] to Ephrin A2
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant full length Human EFNA2 protein (NP_001396) produced in HEK293T cell.
常规说明	Dilute in PBS (pH 7.3) before use.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: PBS, 1% BSA, 50% Glycerol
纯度	Protein A purified
纯化说明	Purified from Mouse ascites fluids by affinity chromatography.
克隆	单克隆
克隆编号	2C6
同种型	IgG2b

应用

Our [Abpromise guarantee](#) covers the use of **ab117941** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/2000. Predicted molecular weight: 24 kDa.

## 靶标


### 序列相似性

Belongs to the ephrin family.

### 细胞定位

Cell membrane.

## 图片



**All lanes :** Anti-Ephrin A2 antibody [2C6] (ab117941) at 1/2000 dilution

**Lane 1 :** HEK293T cells were transfected with the pCMV6-ENTRY control cDNA for 48 hrs and lysed.

**Lane 2 :** HEK293T cells were transfected with the pCMV6-ENTRY EFNA2 cDNA for 48 hrs and lysed.

Lysates/proteins at 5 µg per lane.

**Predicted band size:** 24 kDa

HEK293T cell lysates were generated from transient transfection of the cDNA clone (RC213728)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors