

Product datasheet

Anti-SCCA-2 antibody [2B9] ab117953

1 图像

概述

产品名称	Anti-SCCA-2抗体[2B9]
描述	小鼠单克隆抗体[2B9] to SCCA-2
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant full length Human SERPINB4 produced in HEK293T cells (NP_002965).
阳性对照	HEK293T cell lysate transfected with pCMV6-ENTRY SERPINB4 cDNA.
常规说明	Dilute in PBS (pH7.3) before use. Stable for 12 months from date of receipt.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 48% PBS, 50% Glycerol, 1% BSA
纯度	Protein G purified
纯化说明	ab117953 is purified from Mouse ascites fluid by affinity chromatography.
克隆	单克隆
克隆编号	2B9
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab117953** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500. Predicted molecular weight: 45 kDa.

## 靶标

<b>功能</b>	May act as a protease inhibitor to modulate the host immune response against tumor cells.
<b>组织特异性</b>	Squamous cells.
<b>序列相似性</b>	Belongs to the serpin family. Ov-serpin subfamily.
<b>细胞定位</b>	Cytoplasm. Seems to also be secreted in plasma by cancerous cells but at a low level.

## 图片



**All lanes :** Anti-SCCA-2 antibody [2B9] (ab117953) at 1/500 dilution

**Lane 1 :** HEK293T cell lysate transfected with pCMV6-ENTRY control cDNA

**Lane 2 :** HEK293T cell lysate transfected with pCMV6-ENTRY SERPINB4 cDNA

Lysates/proteins at 5 µg per lane.

**Predicted band size:** 45 kDa

HEK293T cell lysates were generated from transient transfection of the cDNA clone (RC200991)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors

