


Product datasheet

Anti-MRPS36 antibody ab118135

2 图像

概述

产品名称	Anti-MRPS36抗体
描述	小鼠单克隆抗体to MRPS36
宿主	Mouse
经测试应用	适用于: WB, ELISA, Sandwich ELISA
种属反应性	与反应: Recombinant fragment 预测可用于: Human 
免疫原	Recombinant full length protein, corresponding to amino acids 1-103 of Human MRPS36 with proprietary tag (NP_150597).

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Constituent: 99% PBS
纯度	Protein A purified
克隆	单克隆
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab118135** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明


应用	Ab 评论	说明
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 11 kDa. This antibody has only been tested in WB against the recombinant protein used as immunogen. We have no data on the detection of endogenous protein.
ELISA		Use at an assay dependent dilution.
Sandwich ELISA		Use a concentration of 10 µg/ml.

靶标

细胞定位

Mitochondrion.


图片



Anti-MRPS36 antibody (ab118135) at 5 µg/ml
+ recombinant immunogen at 0.2 µg

Predicted band size: 11 kDa

Western blot against tagged recombinant immunogen using ab118135 antibody at 5µg/ml dilution. Predicted band size of immunogen is 38.5kDa.



Detection limit for recombinant tagged MRPS36 is 0.1 ng/ml as a capture antibody.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors