

Product datasheet

Anti-HYAL2 antibody ab118593

2 图像

概述

产品名称	Anti-HYAL2抗体
描述	山羊多克隆抗体to HYAL2
宿主	Goat
经测试应用	适用于: IHC-P
种属反应性	与反应: Human
免疫原	Synthetic peptide: RAGPGPTVTLA-C , corresponding to amino acids 2-12 of Human HYAL2 (NP_003764.3). Run BLAST with Run BLAST with
阳性对照	Human kidney tissue

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 99% Tris buffered saline, 0.5% BSA
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab118593** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		Use a concentration of 3.75 - 7.5 µg/ml.

靶标


功能	Hydrolyzes high molecular weight hyaluronic acid to produce an intermediate-sized product which is further hydrolyzed by sperm hyaluronidase to give small oligosaccharides. Displays very low levels of activity. Associates with and negatively regulates MST1R.
组织特异性	Widely expressed. No expression detected in adult brain.
序列相似性	Belongs to the glycosyl hydrolase 56 family. Contains 1 EGF-like domain.
细胞定位	Cell membrane.

图片



ab118593, at 7.5 µg/ml, staining HYAL2 in formalin fixed, paraffin embedded Human kidney tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-HYAL2 antibody (ab118593)



ab118593, at 7.5 µg/ml, staining HYAL2 in formalin fixed, paraffin embedded Human kidney tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-HYAL2 antibody (ab118593)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors