

Product datasheet

Anti-APBB3 antibody [1G5] ab119017

2 图像

概述

产品名称	Anti-APBB3抗体[1G5]
描述	小鼠单克隆抗体[1G5] to APBB3
宿主	Mouse
经测试应用	适用于: WB, Flow Cyt
种属反应性	与反应: Human
免疫原	Recombinant full length Human APBB3 (NP_573419) produced in HEK293T cells.
阳性对照	HEK293T cell lysate transfected with pCMV6-ENTRY APBB3 cDNA; HEK293T cells transfected with an APBB3 overexpress plasmid

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 48% PBS, 1% BSA, 50% Glycerol
纯度	Protein G purified
纯化说明	ab119017 was purified from Mouse ascites fluids by affinity chromatography.
克隆	单克隆
克隆编号	1G5
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab119017** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/2000. Predicted molecular weight: 53 kDa.

应用	Ab评论	说明
----	------	----

Flow Cyt

1/100.

[ab170190](#) - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

靶标


功能 May modulate the internalization of beta-amyloid precursor protein.

组织特异性 Expressed in various tissues.

序列相似性 Contains 2 PID domains.

Contains 1 WW domain.

图片



All lanes : Anti-APBB3 antibody [1G5]

(ab119017) at 1/2000 dilution


Lane 1 : HEK293T cell lysate transfected with pCMV6-ENTRY control

Lane 2 : HEK293T cell lysate transfected with pCMV6-ENTRY APBB3 cDNA

Lysates/proteins at 5 µg per lane.

Predicted band size: 53 kDa

HEK293T cell lysates were generated from transient transfection of the cDNA clone (RC219752)



HEK293T cells transfected with either an APPB3 overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained using ab119017 at 1/100 dilution, and then analyzed by Flow Cytometry.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors