

Product datasheet

Anti-FAM119A antibody [2D6] ab119186

1 图像

概述

|       |   |
|-------|---|
| 产品名称  | Anti-FAM119A抗体[2D6]   |
| 描述    | 小鼠单克隆抗体[2D6] to FAM119A   |
| 经测试应用 | 适用于: WB   |
| 种属反应性 | 与反应: Human  |
| 免疫原   | Recombinant full length Human FAM119A (NP_660323) produced in HEK293T cell. |
| 阳性对照  | HEK293T cells transfected with the pCMV6-ENTRY FAM119A                      |

性能

|      |  |
|------|--|
| 形式   | Liquid   |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | pH: 7.30<br>Preservative: 0.02% Sodium azide<br>Constituents: 48% PBS, 1% BSA, 50% Glycerol    |
| 纯度   | IgG fraction   |
| 纯化说明 | Purified from mouse ascites fluids by affinity chromatography                                  |
| 克隆   | 单克隆  |
| 克隆编号 | 2D6  |
| 同种型  | IgG1   |

应用

Our [Abpromise guarantee](#) covers the use of **ab119186** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明  |
|----|------|---|
| WB |      | 1/2000. Predicted molecular weight: 25 kDa. |

图标


## 序列相似性

Belongs to the FAM119 family.

## 细胞定位

Membrane.

## 图片



**All lanes :** Anti-FAM119A antibody [2D6] (ab119186) at 1/2000 dilution

**Lane 1 :** HEK293T cells transfected with the pCMV6-ENTRY control

**Lane 2 :** HEK293T cells transfected with the pCMV6-ENTRY FAM119A

Lysates/proteins at 5 µg per lane.

**Predicted band size :** 25 kDa

HEK293T cell lysates were generated from transient transfection of the cDNA clone (RC202889)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors