

Product datasheet

Anti-TRAF5 antibody [55A219] ab12123

1 References 1 图像

概述

产品名称	Anti-TRAF5抗体[55A219]
描述	小鼠单克隆抗体[55A219] to TRAF5
宿主	Mouse
特异性	ab12123 reacts with human TRAF5.
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Fusion protein corresponding to amino acid residue(s) 77 to 186 of human TRAF5.
阳性对照	HeLa whole-cell extract.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: PBS
纯度	Protein G purified
纯化说明	Affinity purified.
克隆	单克隆
克隆编号	55A219
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab12123** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 - 2 µg/ml. Detects a band of approximately 66 kDa (predicted molecular weight: 66.8 kDa).

靶标

功能	Adapter protein and signal transducer that links members of the tumor necrosis factor receptor family to different signaling pathways by association with the receptor cytoplasmic domain and kinases. Mediates activation of NF-kappa-B and probably JNK. Seems to be involved in apoptosis.
组织特异性	Expressed in spleen, thymus, prostate, testis, ovary, small intestine, colon, and peripheral blood.
序列相似性	Belongs to the TNF receptor-associated factor family. A subfamily. Contains 1 MATH domain. Contains 1 RING-type zinc finger. Contains 2 TRAF-type zinc fingers.
结构域	The MATH/TRAF domain binds to receptor cytoplasmic domains.
细胞定位	Cytoplasm.

图片



TRAF5 detection by Western blot. The analysis of 10 µg of HeLa cell extract by using ab12123 at a 2 µg/ml dilution.

TRAF5 detection by Western blot. The analysis of 10 µg of HeLa cell extract by using ab12123 at a 2 µg/ml dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors