

Product datasheet

Anti-BANF2 antibody ab122950

1 图像

概述

产品名称	Anti-BANF2抗体
描述	兔多克隆抗体to BANF2
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse 预测可用于: Rat, Rabbit, Horse, Guinea pig, Cow, Cat, Dog, Human, Pig
免疫原	Synthetic peptide corresponding to a region within N terminal amino acids 2-51 (DDMSPRLRAF LSEPIGEKDV AWVDGISREL AINLVTKGFN KAYVLLGQFL) of Mouse BANF2 (NP_997158). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Mouse small intestine lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituents: 97% PBS, 2% Sucrose
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab122950** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 10 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

功能	May play a role in BANF1 regulation and influence tissue-specific roles of BANF1.
组织特异性	Expressed strongly in testis and pancreas. Also detected in brain, colon, liver, lung, ovary, placenta, prostate, small intestine, spleen and thymus. Not detected in heart, kidney and skeletal muscle.
序列相似性	Belongs to the BAF family.
细胞定位	Nucleus. Cytoplasm.

图片



Anti-BANF2 antibody (ab122950) at 1 μ g/ml +
Mouse small intestine lysate at 10 μ g

Predicted band size: 10 kDa

Gel concentration: 10-20%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors