

Product datasheet

Anti-PAP2D antibody ab124177

1 图像

概述

产品名称	Anti-PAP2D抗体
描述	兔多克隆抗体to PAP2D
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region within internal sequence amino acids 77-106 of Human PAP2D (NP_001010861.1, NP_001032394.1).
阳性对照	HepG2 cell line lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium azide Constituent: 99% PBS
纯度	Immunogen affinity purified
纯化说明	ab124177 is purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab124177** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/100 - 1/500. Predicted molecular weight: 35 kDa.

标签


---

组织特异性	Isoform 1 is expressed in brain, lung, kidney and colon. Isoform 2 is expressed in placenta, skeletal muscle and kidney.
序列相似性	Belongs to the PA-phosphatase related phosphoesterase family.
细胞定位	Membrane.

---

## 图片

---



Anti-PAP2D antibody (ab124177) at 1/100 dilution + HepG2 cell line lysate at 35 µg

**Predicted band size:** 35 kDa

Western blot - Anti-PAP2D antibody (ab124177)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors