

Product datasheet

Anti-LY108 antibody ab124255

1 图像

概述

产品名称	Anti-LY108抗体
描述	兔多克隆抗体to LY108
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region within C-terminal amino acids 276-305 of Human LY108 (NP_443163.1, NP_001171644.1).
阳性对照	T47D cell line lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium azide Constituent: 99% PBS
纯度	Immunogen affinity purified
纯化说明	ab124255 is purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab124255** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		1/100 - 1/500. Predicted molecular weight: 37 kDa.

靶标

功能	Triggers cytolytic activity only in natural killer cells (NK) expressing high surface densities of natural cytotoxicity receptors.
----	--

组织特异性	Expressed by all (resting and activated) natural killer cells (NK), T- and B-lymphocytes.
序列相似性	Contains 1 Ig-like (immunoglobulin-like) domain.
细胞定位	Cell membrane.

Anti-LY108 antibody 图像



Anti-LY108 antibody (ab124255) at 1/100 dilution + T47D cell line lysate at 35 μ g

Predicted band size : 37 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors