


Product datasheet

Anti-ADRM1 antibody ab124645

2 图像

概述

产品名称	Anti-ADRM1抗体
描述	兔多克隆抗体to ADRM1
宿主	Rabbit
经测试应用	适用于: IHC-P
种属反应性	与反应: Mouse, Human 预测可用于: Chimpanzee, Rhesus monkey, Orangutan 
免疫原	Synthetic peptide corresponding to a region between residues 357 and 407 of Human ADRM1 (NP_008933.2).
阳性对照	Human osteosarcoma, ovarian carcinoma and prostate carcinoma tissues; Mouse teratoma tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium azide Constituents: 99% Tris buffered saline, 0.1% BSA
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab124645** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		1/100 - 1/500. Epitope exposure is recommended. Epitope exposure with citrate buffer will enhance staining.

靶标


功能	Functions as a proteasomal ubiquitin receptor. Recruits the deubiquitinating enzyme UCHL5 at the 26S proteasome and promotes its activity.
序列相似性	Belongs to the ADRM1 family. Contains 1 PH domain.
结构域	The PH domain mediates interactions with PSMD1 and ubiquitin. Preferential binding to the proximal subunit of K48-linked diubiquitin allows UCHL5 access to the distal subunit.
细胞定位	Cytoplasm. Nucleus.

图片



ab124645, at 1/250 dilution, staining ADRM1 in formalin-fixed, paraffin-embedded Human prostate carcinoma by Immunohistochemistry with DAB staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ADRM1 antibody (ab124645)



ab124645, at 1/250 dilution, staining ADRM1 in formalin-fixed, paraffin-embedded Mouse teratoma tissue by Immunohistochemistry with DAB staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ADRM1 antibody (ab124645)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors