

Product datasheet

Anti-RRM1 antibody ab125713

2 图像

概述

产品名称	Anti-RRM1 抗体
描述	兔多克隆抗体 to RRM1
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human
免疫原	Recombinant fragment containing a sequence corresponding to a region within amino acids 246-765 of Human RRM1.(Uniprot:P23921).
阳性对照	293T cell lysate, A431 cell lysate, Jurkat cell lysate and Human Colon cancer tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 78% PBS, 20% Glycerol, 1% BSA
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab125713** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 90 kDa.

应用	Ab评论	说明
IHC-P		1/100 - 1/1000. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## 靶标


功能	Provides the precursors necessary for DNA synthesis. Catalyzes the biosynthesis of deoxyribonucleotides from the corresponding ribonucleotides.
通路	Genetic information processing; DNA replication.
序列相似性	Belongs to the ribonucleoside diphosphate reductase large chain family. Contains 1 ATP-cone domain.
细胞定位	Cytoplasm.

## 图片



Anti-RRM1 antibody (ab125713) at 1/2000 dilution + 293T lysate at 30  $\mu$ g

**Predicted band size:** 90 kDa



ab125713, at 1/500 dilution, staining RRM1 in Formalin-fixed, Paraffin-embedded Human colon cancer tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-RRM1 antibody (ab125713)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors