

Product datasheet

Anti-Bacillus anthracis protective antigen antibody ab13809

1 References 1 图像

概述

产品名称	Anti-Bacillus anthracis protective antigen抗体
描述	兔多克隆抗体 to Bacillus anthracis protective antigen
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	Anthrax PA produced by Bacillus anthracis
免疫原	Synthetic peptide corresponding to amino acids 186-192 of Anthrax PA

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.05% Sodium Azide Constituents: 0.2% Gelatin, PBS
纯度	Protein G purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab13809** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/1000. Detects a band of approximately 83 kDa.

靶标


相关性 Anthrax toxins are composed of three distinct proteins, a protective antigen (PA), a lethal factor

(LF) and an edema factor (EF). None of these proteins are toxic by themselves and several studies indicate that the anthrax toxin has the familiar A-B enzymatic binding structure, with PA acting as the binding domain and EF and/or LF acting as the active fragments.

细胞定位

Secreted

图片



Western blot analysis of Anthrax PA in recombinant protein using ab13809 (lane B) at 1:1000 dilution. Lane A shows the pre-bleed.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors