

Product datasheet

Anti-Endostatin antibody [1837-46] ab15685

1 图像

概述

<b>产品名称</b>	Anti-Endostatin抗体[1837-46]
<b>描述</b>	小鼠单克隆抗体[1837-46] to Endostatin
<b>宿主</b>	Mouse
<b>经测试应用</b>	<b>适用于:</b> IHC-Fr, IHC-P
<b>种属反应性</b>	<b>与反应:</b> Human
<b>免疫原</b>	Recombinant fragment (Human) (C terminal).
<b>常规说明</b>	Clone 1837-46 is directed against the endostatin portion (the carboxy-terminal domain) of collagen XVIII. Endostatin is a powerful inhibitor of endothelial cell proliferation and angiogenesis and shows widespread distribution, mainly in vessel walls and some other basement membrane zones.

性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
<b>存储溶液</b>	Preservative: 0.1% Sodium Azide Constituents: PBS, pH 7.2
<b>纯度</b>	Immunogen affinity purified
<b>纯化说明</b>	Purified from tissue culture supernatant.
<b>Primary antibody说明</b>	Clone 1837-46 is directed against the endostatin portion (the carboxy-terminal domain) of collagen XVIII. Endostatin is a powerful inhibitor of endothelial cell proliferation and angiogenesis and shows widespread distribution, mainly in vessel walls and some other basement membrane zones.
<b>克隆</b>	单克隆
<b>克隆编号</b>	1837-46
<b>同种型</b>	IgG2b
<b>轻链类型</b>	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab15685** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-Fr		Use a concentration of 10 µg/ml. Fixation: acetone 10 minutes.
IHC-P		Use a concentration of 10 µg/ml. Microwave irradiation (2 x 7min) in citrate buffer is recommended for antigen retrieval.

## 靶标

<b>功能</b>	COLA18A probably plays a major role in determining the retinal structure as well as in the closure of the neural tube. Endostatin potently inhibits endothelial cell proliferation and angiogenesis. May inhibit angiogenesis by binding to the heparan sulfate proteoglycans involved in growth factor signaling.
<b>组织特异性</b>	Present in multiple organs with highest levels in liver, lung and kidney.
<b>疾病相关</b>	Defects in COL18A1 are a cause of Knobloch syndrome (KNO) [MIM:267750]. KNO is an autosomal recessive disorder defined by the occurrence of high myopia, vitreoretinal degeneration with retinal detachment, macular abnormalities and occipital encephalocele.
<b>序列相似性</b>	Belongs to the multiplexin collagen family. Contains 1 FZ (frizzled) domain. Contains 1 TSP N-terminal (TSPN) domain.
<b>翻译后修饰</b>	Prolines at the third position of the tripeptide repeating unit (G-X-Y) are hydroxylated in some or all of the chains.
<b>细胞定位</b>	Secreted > extracellular space > extracellular matrix.

## 图片



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Endostatin antibody [1837-46] (ab15685)

IHC image of ab15685 staining in human colon formalin fixed paraffin embedded tissue section, performed on a Leica Bond™ system using the standard protocol F. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH6, epitope retrieval solution 1) for 20 mins. The section was then incubated with ab15685, 1 µg/ml, for 15 mins at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

For other IHC staining systems (automated and non-automated) customers should optimize variable parameters such as antigen retrieval conditions, primary antibody concentration and antibody incubation times.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors