

Product datasheet

Anti-CD37 antibody [IPO-24] ab19701

概述

产品名称	Anti-CD37抗体[IPO-24]
描述	小鼠单克隆抗体[IPO-24] to CD37
宿主	Mouse
经测试应用	适用于: IP, ELISA, IHC-P, IHC-Fr
种属反应性	与反应: Human
免疫原	Tissue/ cell preparation (B cell line).

性能

形式	Liquid
存放说明	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Mix suspension via gentle inversion before use.
存储溶液	Preservative: None Constituents: PBS, pH 7.2
纯度	Protein A purified
克隆	单克隆
克隆编号	IPO-24
同种型	IgG2b

应用

Our [Abpromise guarantee](#) covers the use of **ab19701** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IP		Use at an assay dependent concentration.
ELISA		1/1000.
IHC-P		Use at an assay dependent concentration.

应用	Ab评论	说明
IHC-Fr		Use at an assay dependent concentration.

靶标

组织特异性	B-lymphocytes.
序列相似性	Belongs to the tetraspanin (TM4SF) family.
细胞定位	Membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors