

Product datasheet

Anti-Transglutaminase 2 antibody [2A1] ab21260

1 References 1 图像

概述

产品名称	Anti-Transglutaminase 2抗体[2A1]
描述	小鼠单克隆抗体[2A1] to Transglutaminase 2
宿主	Mouse
经测试应用	适用于: ELISA, WB, IP
种属反应性	与反应: Human 不与反应: Mouse, Rat
免疫原	Recombinant human protein purified from E.coli.
阳性对照	U87MG cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 50% Glycerol, PBS
纯度	Protein G purified
克隆	单克隆
克隆编号	2A1
同种型	IgG1
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab21260** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		

应用	Ab评论	说明
----	------	----

WB

IP

#### 应用说明

ELISA: Use at an assay dependent dilution.  
 IP: Use at a concentration of 1 - 2 µg/ml.  
 WB: 1/2000. Detects a band of approximately 78 kDa.


Not yet tested in other applications.  
 Optimal dilutions/concentrations should be determined by the end user.

#### 靶标

**功能** Catalyzes the cross-linking of proteins and the conjugation of polyamines to proteins.

**序列相似性** Belongs to the transglutaminase superfamily. Transglutaminase family.

#### 图片



Anti-Transglutaminase 2 antibody [2A1] (ab21260) at 1/2000 dilution + U87MG cell lysate

**Observed band size:** 78 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors