

## Product datasheet

# Anti-TRIM29 antibody ab22207

★★★★☆ 1 Abreviews 2 References 1 图像

### 概述

产品名称	Anti-TRIM29抗体
描述	兔多克隆抗体to TRIM29
宿主	Rabbit
特异性	ab22207 reacts with TRIM29.
经测试应用	<b>适用于:</b> IHC-P, ELISA, WB
种属反应性	<b>与反应:</b> Human <b>预测可用于:</b> Mouse, Rat, Dog, Chimpanzee
免疫原	Synthetic peptide: SHYRPFYVNK GNGIG , corresponding to C terminal amino acids 569-583 of Human TRIM29. <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
阳性对照	Testis tissue

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.05% Sodium azide Constituents: PBS, 0.2% Gelatin
纯度	Protein G purified
克隆	多克隆
同种型	IgG

### 应用


Our [Abpromise guarantee](#) covers the use of **ab22207** in the following tested applications.  
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		Use at an assay dependent concentration.
ELISA		1/100 - 1/1000.
WB	★★★★☆	Use at an assay dependent concentration.

## 靶标

功能	It is able to complement the radiosensitivity defect of an ataxia telangiectasia (AT) fibroblast cell line.
组织特异性	Expressed in placenta, prostate and thymus.
序列相似性	Contains 1 B box-type zinc finger.
翻译后修饰	Constitutively phosphorylated by PKC on serine/threonine in A431 cells.
细胞定位	Cytoplasm. Colocalizes with intermediate filaments.

## 图片



Immunohistochemical staining with ab22207 in normal testis tissue at 1:100 dilution.

Immunohistochemistry (Paraffin-embedded sections)  
- TRIM29 antibody (ab22207)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors