

Product datasheet

Anti-TRIM29 antibody ab22207

★★★★☆ 1 Abreviews 2 References 1 图像

概述

产品名称	Anti-TRIM29抗体
描述	兔多克隆抗体to TRIM29
宿主	Rabbit
特异性	ab22207 reacts with TRIM29.
经测试应用	适用于: IHC-P, ELISA, WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Dog, Chimpanzee
免疫原	Synthetic peptide: SHYRPFYVNK GNGIG , corresponding to C terminal amino acids 569-583 of Human TRIM29. Run BLAST with Run BLAST with
阳性对照	Testis tissue

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.05% Sodium azide Constituents: PBS, 0.2% Gelatin
纯度	Protein G purified
克隆	多克隆
同种型	IgG

应用


Our [Abpromise guarantee](#) covers the use of **ab22207** in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		Use at an assay dependent concentration.
ELISA		1/100 - 1/1000.
WB	★★★★☆	Use at an assay dependent concentration.

靶标

功能	It is able to complement the radiosensitivity defect of an ataxia telangiectasia (AT) fibroblast cell line.
组织特异性	Expressed in placenta, prostate and thymus.
序列相似性	Contains 1 B box-type zinc finger.
翻译后修饰	Constitutively phosphorylated by PKC on serine/threonine in A431 cells.
细胞定位	Cytoplasm. Colocalizes with intermediate filaments.

图片



Immunohistochemical staining with ab22207 in normal testis tissue at 1:100 dilution.

Immunohistochemistry (Paraffin-embedded sections)
- TRIM29 antibody (ab22207)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors