


# Anti-ACTH antibody ab2407

## 概述

产品名称	Anti-ACTH抗体
描述	兔多克隆抗体to ACTH
宿主	Rabbit
特异性	This antibody reacts with human ACTH. This antibody cross reacts with ACTH in a wide variety of mammalian species. Traces of contaminating antibodies have been removed by solid phase absorption.
经测试应用	适用于: IHC-P, IHC-Fr
种属反应性	与反应: Human 预测可用于: a wide range of other species 
免疫原	Synthetic human ACTH 1-24.
阳性对照	Pituitary

## 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	unknown
轻链类型	unknown

## 应用

Our [Abpromise guarantee](#) covers the use of **ab2407** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		
IHC-Fr		
<b>应用说明</b>	This antibody may be diluted to a titer of 1:50-1:250 in an ABC method. We suggest an incubation period of 30 minutes at room temperature. Suitable for formalin/paraffin tissues.	
<b>靶标</b>		
<b>相关性</b>	ACTH occurs in cells of the anterior pituitary and in neurons in brain. It regulates the corticosteroid production in the adrenal cortex. Beta endorphin and Met enkephalin are endogenous opiates. MSH (melanocyte stimulating hormone) increases the pigmentation of skin by increasing melanin production in melanocytes.	
<b>细胞定位</b>	Secreted	

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors