

Product datasheet

Anti-Qa 2 antibody [695H1-9-9] ab25366

概述

产品名称	Anti-Qa 2抗体[695H1-9-9]
描述	小鼠单克隆抗体[695H1-9-9] to Qa 2
宿主	Mouse
特异性	This antibody reacts with a Qa2 determined alloantigen.
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse
免疫原	Tissue/ cell preparation (Mouse): C3H.SW skin graft and splenocytes.
常规说明	

This antibody may be useful for Flow Cytometry, studies on complement mediated cytotoxicity and modulation of Qa2 expression from peripheral T lymphocytes.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 8.20 Constituent: 100% Borate buffered saline
纯度	IgG fraction
Primary antibody说明	This antibody may be useful for Flow Cytometry, studies on complement mediated cytotoxicity and modulation of Qa2 expression from peripheral T lymphocytes.
克隆	单克隆
克隆编号	695H1-9-9
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab25366** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Flow Cyt		Use 1µg for 10 ⁶ cells. ab170191 - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody.

靶标

相关性 Qa antigens are 'non classical' MHC class I molecules of mice. Qa2 is present on restricted functional lymphocyte subpopulations. The protein product of the preimplantation embryo development (Ped) gene, Qa2 influences the rate of embryonic cleavage during the preimplantation stages of development.

细胞定位 Plasma membrane

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors