

Product datasheet

Anti-TBC1D10B antibody ab28632

1 图像

概述

产品名称	Anti-TBC1D10B抗体
描述	兔多克隆抗体to TBC1D10B
宿主	Rabbit
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Human
免疫原	A region within a synthetic peptide: EKEREKQEKE RQKQEKKAAQG RKLSLRRKAD GPPGPHDGGD RPSAEARQDA, corresponding to C terminal amino acids 757-806 of Human TBC1D10B <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
阳性对照	HepG2 cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab28632** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
ELISA		1/12500.

应用	Ab评论	说明
WB		Use a concentration of 1.25 µg/ml. Predicted molecular weight: 61 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

**相关性** TBC1D10B, which belongs to the TBC1 domain family, member 10B, was experimentally validated as a HeLa cell nuclei protein, currently with unknown function . TBC/rabGAP domain-containing microvillar protein binds to first PDZ domains of EBP50 and E3KARP and EBP50-PDZ and may be potentially involved in regulating membrane traffic.

## 图片



Anti-TBC1D10B antibody (ab28632) at 1.25 µg/ml + HepG2 Cell lysate 10ug.

**Predicted band size:** 61 kDa

**Observed band size:** 59 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors