

Product datasheet

Anti-HSV2 gG Envelope antibody [M680910] ab32804

概述

产品名称	Anti-HSV2 gG Envelope抗体[M680910]
描述	小鼠单克隆抗体[M680910] to HSV2 gG Envelope
宿主	Mouse
特异性	Recognizes glycoprotein G (gG) of Herpes Simplex Virus Type 2 only. Reactive in IFA of infected cells, Western blot and immunoprecipitation of specific viral protein. The clone number has been updated from (8.F.141) to (M680910) both clone numbers name the same antibody clone.
经测试应用	适用于: WB, ICC/IF, IP
种属反应性	Herpes Simplex Virus II
免疫原	Herpes Simplex Virus II Glycoprotein G (gG).

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Add glycerol to a final volume of 50% for extra stability and aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.01% Sodium Azide Constituents: PBS, pH 7.4
纯度	IgG fraction
克隆	单克隆
克隆编号	M680910
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab32804** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use at an assay dependent dilution. Predicted molecular weight: 8 kDa.

应用	Ab评论	说明
ICC/IF		Use at an assay dependent dilution.
IP		Use at an assay dependent dilution.

靶标

相关性 HSV2 is a member of the herpes virus family, Herpesviridae. All herpes viruses are morphologically identical. They have a large double stranded DNA genome. The virion consists of an icosahedral nucleocapsid which is surrounded by a lipid bilayer envelope. Between the capsid and the envelope is an amorphous layer of proteins, termed the tegument. Glycoprotein G (gG) is one of seven envelope glycoproteins in HSV2.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors