

Product datasheet

Anti-AU1 tag antibody (HRP) ab3402

概述

| | |
|-------|--|
| 产品名称 | Anti-AU1 tag抗体(HRP) |
| 描述 | 山羊多克隆抗体to AU1 tag (HRP) |
| 宿主 | Goat |
| 偶联物 | HRP |
| 经测试应用 | 适用于: ELISA, ICC, WB |
| 免疫原 | Full length native AU1 (purified): DTYRYI conjugated to KLH. |

[Run BLAST with](#)

[Run BLAST with](#)

常规说明 Molar enzyme/antibody protein ratio is 4:1.

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C. |
| 存储溶液 | pH: 7.00 Constituents: 1.19% HEPES, 0.58% Sodium chloride, 0.2% BSA |
| 纯度 | Immunogen affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab3402** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|-------|------|---------------------|
| ELISA | | 1/10000 - 1/100000. |
| ICC | | 1/200 - 1/500. |

| 应用 | Ab评论 | 说明 |
|----|------|--------------------|
| WB | | 1/10000 - 1/30000. |

靶标

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors