

Product datasheet

Anti-HSV2 gC antibody [M690403] ab46860

概述

产品名称	Anti-HSV2 gC抗体[M690403]
描述	小鼠单克隆抗体[M690403] to HSV2 gC
特异性	Reactive with Glycoprotein C (gC-2) of Herpes Simplex Virus Type 2 only.
经测试应用	适用于: IP, ICC/IF
种属反应性	Reacts with Herpes Simplex Virus 2. Not yet tested in other species.
免疫原	Herpes Simplex Virus 2 Glycoprotein C (gC-2).

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS, pH 7.4
克隆	单克隆
克隆编号	M690403
同种型	IgG2a

应用

Our [Abpromise guarantee](#) covers the use of **ab46860** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IP		Use at an assay dependent dilution.
ICC/IF		Use at an assay dependent dilution.

靶标

**相关性** All herpes viruses are morphologically identical: they have a large double stranded DNA genome. The virion consists of an icosahedral nucleocapsid which is surrounded by a lipid

bilayer envelope. Between the capsid and the envelope is an amorphous layer of proteins, termed the tegument. A characteristic of all herpes viruses is that, following primary infection, the virus establishes a latent infection in the host and may be reactivate at any stage. Reactivation is frequently, but not always, associated with further disease. HSV2 is usually associated with genital lesions. There are seven external glycoproteins in HSV1 and 2: gH, gB, gC, gG, gD, gI, and gE.

**细胞定位**

External viral glycoprotein

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

**Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors