

Product datasheet

Anti-GPCR GPR44 antibody ab46931

1 图像

概述

<b>产品名称</b>	Anti-GPCR GPR44抗体
<b>描述</b>	兔多克隆抗体to GPCR GPR44
<b>经测试应用</b>	<b>适用于:</b> WB
<b>种属反应性</b>	<b>与反应:</b> Mouse, Rat, Human <b>预测可用于:</b> Chimpanzee, Macaque monkey 
<b>免疫原</b>	Synthetic peptide: DDSELGGAGSSRRRRTSSTA , corresponding to C terminal amino acids 329-348 of Human GPCR GPR44. <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
<b>阳性对照</b>	GPCR GPR44-transfected Jurkat cells.

性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>存储溶液</b>	Preservative: 0.02% Sodium Azide Constituents: PBS
<b>纯度</b>	Protein G purified
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab46931** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 66 kDa (predicted molecular weight: 43 kDa).

## 靶标


---

功能	Orphan receptor.
序列相似性	Belongs to the G-protein coupled receptor 1 family.
细胞定位	Cell membrane.

---

## 图片

---



Anti-GPCR GPR44 antibody (ab46931) at  
1/500 dilution + GPR44 transfected Jurkat  
cells

**Predicted band size** : 43 kDa

**Observed band size** : 66 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors