


Product datasheet

Anti-Cortactin antibody ab47582

1 图像

概述

产品名称	Anti-Cortactin抗体
描述	兔多克隆抗体to Cortactin
宿主	Rabbit
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	A synthesized non-phosphopeptide derived from human Cortactin around the phosphorylation site of tyrosine 421 (P-IY <sub>p</sub> -E-D)
阳性对照	Extracts from HeLa cells

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS, 150mM Sodium chloride, pH 7.4
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab47582** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 85 kDa (predicted molecular weight: 62 kDa).
ELISA		1/20000.

## 靶标

<b>功能</b>	Contributes to the organization of the actin cytoskeleton and cell structure. Plays a role in the regulation of cell migration. Plays a role in the invasiveness of cancer cells, and the formation of metastases.
<b>序列相似性</b>	Contains 7 cortactin repeats. Contains 1 SH3 domain.
<b>结构域</b>	The SH3 motif may mediate binding to the cytoskeleton.
<b>翻译后修饰</b>	Tyrosine phosphorylation in transformed cells may contribute to cellular growth regulation and transformation.
<b>细胞定位</b>	Cytoplasm > cytoskeleton. Cell projection > lamellipodium. Cell projection > ruffle. Associated with membrane ruffles and lamellipodia.

## 图片



**All lanes :** Anti-Cortactin antibody (ab47582)  
at 1/500 dilution

**Lane 1 :** Extracts from HeLa cells treated with hydrogen peroxide.

**Lane 2 :** Extracts from HeLa cells treated with hydrogen peroxide and immunogenic peptide.

Lysates/proteins at 30 µg per lane.

**Predicted band size:** 62 kDa

**Observed band size:** 85 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- 
- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors