

## Product datasheet

# Anti-Endothelin 2 antibody ab47923

### 概述

产品名称	Anti-Endothelin 2抗体
描述	兔多克隆抗体to Endothelin 2
特异性	This antibody can cross react with Endothelin-1 (20%) and Endothelin-3 (less than 10%).
经测试应用	适用于: RIA, IP, ELISA
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Cow, Cat, Dog, Ferret
免疫原	Synthetic peptide: CSCSSWLDKECVYFCHLDIIW conjugated to KLH via Cys, corresponding to amino acids 49-69 of Human Endothelin 2. <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 50% Glycerol, PBS, pH 7.5
纯度	Protein G purified
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab47923** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
RIA		Use at an assay dependent dilution.
IP		Use at an assay dependent dilution.

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.

## 靶标

功能	Endothelins are endothelium-derived vasoconstrictor peptides.
组织特异性	Expressed in lung, but not in placental stem villi vessels or cultured placental villi smooth muscle cells.
序列相似性	Belongs to the endothelin/sarafotoxin family.
细胞定位	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors