

Product datasheet

Anti-HOXB7 antibody ab48829

1 图像

概述

产品名称	Anti-HOXB7抗体
描述	兔多克隆抗体to HOXB7
宿主	Rabbit
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Mouse 预测可用于: Rat, Cow, Dog, Human, Chimpanzee 
免疫原	A region within synthetic peptide: GAKEQRDSDL AAESNFRYP WMRSSGPDRK RGRQTYTRYQ TLELEKEFHY, corresponding to amino acids 109-158 of Human HOXB7 Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	<div style="border: 1px solid #ccc; padding: 5px; display: inline-block;"> <p>购买相配的WB阳性对照 Recombinant Human HOXB7 protein ></p> </div> large intestine lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab48829** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
ELISA		1/312500.
WB		Use a concentration of 1.25 µg/ml. Detects a band of approximately 24 kDa (predicted molecular weight: 24 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

功能	Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.
序列相似性	Belongs to the Antp homeobox family. Contains 1 homeobox DNA-binding domain.
细胞定位	Nucleus.

图片



Anti-HOXB7 antibody (ab48829) at 1.25 µg/ml
+ large intestine lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size: 24 kDa

Observed band size: 24 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors