

Product datasheet

Anti-Annexin IV antibody [NO.73] ab50142

概述

产品名称	Anti-Annexin IV抗体[NO.73]
描述	小鼠单克隆抗体[NO.73] to Annexin IV
宿主	Mouse
经测试应用	适用于: Sandwich ELISA, ELISA
种属反应性	与反应: Human
免疫原	Full length native protein (purified)

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: 0.15M Sodium chloride, 0.15mM PBS, pH 7.2
纯度	Protein A purified
克隆	单克隆
克隆编号	NO.73
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab50142** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Sandwich ELISA		
ELISA		

应用说明	ELISA: Use at a concentration of 1 - 10 µg/ml. sELISA: Can be paired for ELISA with Mouse monoclonal [NO.50] to Annexin IV (ab50141) . ab50142 can be used as a detection antibody in a human Annexin IV ELISA combination with
------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

the monoclonal capture antibody [ab50141](#). Use at an assay dependent dilution.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

功能	Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.
序列相似性	Belongs to the annexin family. Contains 4 annexin repeats.
结构域	A pair of annexin repeats may form one binding site for calcium and phospholipid.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors