

## Product datasheet

# Anti-ZNF253 antibody ab50992

### 1 图像

#### 概述

产品名称	Anti-ZNF253抗体
描述	兔多克隆抗体to ZNF253
宿主	Rabbit
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human
免疫原	A region within synthetic peptide: FTASSTLTTH KRIHTGEKPYKCEECGKAFN WSSDLNKHKK IHIERKPYIV, corresponding to C terminal amino acids 433-482 of Human ZNF253 <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
阳性对照	HepG2 cell lysate

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab50992** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1.25 µg/ml. Predicted molecular weight: 32 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

应用	Ab评论	说明
----	------	----


ELISA Use at an assay dependent concentration.  
ELISA titre using peptide based assay: 1:312500.

## 靶标

**相关性** ZNF253 may function as a transcription factor, possibly by exerting a transcriptional repression activity.

**细胞定位** Nuclear

## 图片



**Lane 1** : MW marker

**Lane 2** : Anti-ZNF253 antibody (ab50992) at 1.25 µg/ml

**Lane 2** : HepG2 cell lysate at 10 µg

### Secondary

**Lane 2** : HRP conjugated anti-Rabbit IgG at 1/50000 dilution

**Predicted band size:** 32 kDa

Gel concentration: 12%

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors