

Product datasheet

Anti-ALDH7A1 antibody ab51029

1 References 3 图像

概述

产品名称	Anti-ALDH7A1抗体
描述	兔多克隆抗体to ALDH7A1
经测试应用	适用于: WB, IP, IHC-P, ICC/IF
种属反应性	与反应: Mouse, Human 不与反应: Rat
免疫原	Synthetic peptide corresponding to residues near the N terminus of human ALDH7A1.
阳性对照	Fetal liver lysate or HepG2 cells.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	PBS 49%,Sodium azide 0.01%,Glycerol 50%,BSA 0.05%
纯度	Protein A purified
克隆	多克隆
同种型	IgG


应用

Our [Abpromise guarantee](#) covers the use of **ab51029** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/10000. Detects a band of approximately 50 kDa (predicted molecular weight: 55 kDa).
IP		1/40.
IHC-P		Use at an assay dependent concentration.
ICC/IF		1/250 - 1/500.

图片




Anti-ALDH7A1 antibody (ab51029) at
1/10000 dilution + Fetal liver lysate at 10 μ g

Secondary

Goat anti-rabbit HRP labeled at 1/2000
dilution

Predicted band size : 55 kDa

Observed band size : 50 kDa



Ab51029 (1:100) staining human ALDH7A1 in
human liver by immunohistochemistry using
paraffin embedded tissue.



Imunofluorescent staining of HEPG2 cells
using ab51029 (1:250).

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors