

Product datasheet

Anti-S100A9 + Calprotectin (S100A8/A9 complex) antibody [MAC387], prediluted ab51934

概述

产品名称	Anti-S100A9 + Calprotectin (S100A8/A9 complex)抗体[MAC387], prediluted
描述	小鼠单克隆抗体[MAC387] to S100A9 + Calprotectin (S100A8/A9 complex), prediluted
宿主	Mouse
特异性	This antibody is specific for neutrophilic granulocytes, monocytes and reactive macrophages
经测试应用	适用于: IHC-P
种属反应性	与反应: Human
免疫原	Tissue/ cell preparation - monocyte preparation
阳性对照	Spleen
常规说明	This antibody is suitable for the detection of tissue specific histiocytes / macrophages in different tissues.

性能

形式	Prediluted
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: Stabilising phosphate buffer, pH 7.3
Primary antibody说明	This antibody is suitable for the detection of tissue specific histiocytes / macrophages in different tissues.
克隆	单克隆
克隆编号	MAC387
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab51934** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		Use at an assay dependent concentration. 100 µl working solution / test. Pre-treatment with pronase 0.1% or unmasking fluid G is recommended.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors