

Product datasheet

Anti-ID2 antibody, prediluted ab52116

1 图像

概述

产品名称	Anti-ID2抗体, prediluted
描述	兔多克隆抗体到 ID2, prediluted
宿主	Rabbit
经测试应用	适用于: IHC-P
种属反应性	与反应: Human
免疫原	A synthetic peptide from the C-terminus of human Id2 protein
阳性对照	Human breast carcinoma

性能

形式	Prediluted
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: 1% BSA, 50mM Tris, pH 7.6
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab52116** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		


应用说明 IHC-P: Ready to use. Perform heat mediated antigen retrieval before commencing with IHC staining protocol.

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

靶标

功能	ID (inhibitor of DNA binding) HLH proteins lack a basic DNA-binding domain but are able to form heterodimers with other HLH proteins, thereby inhibiting DNA binding. ID-2 may be an inhibitor of tissue-specific gene expression.
组织特异性	Highly expressed in early fetal tissues, including those of the central nervous system.
序列相似性	Contains 1 basic helix-loop-helix (bHLH) domain.
发展阶段	Found in most early fetal tissues but not in the corresponding mature tissues.
细胞定位	Cytoplasm. Nucleus.

图片



ab52116, staining human ID2 in breast carcinoma tissue by Immunohistochemistry, Formalin fixed Paraffin embedded tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - ID2 antibody, prediluted (ab52116)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors