

Product datasheet

Anti-Dlx2 antibody ab54673

2 图像

概述

| | |
|-------|--|
| 产品名称 | Anti-Dlx2抗体 |
| 描述 | 小鼠单克隆抗体to Dlx2 |
| 宿主 | Mouse |
| 经测试应用 | 适用于: WB, IHC-P |
| 种属反应性 | 与反应: Mouse, Human |
| 免疫原 | Recombinant fragment, corresponding to amino acids 1-110 of Human Dlx2 |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: None PBS, pH 7.2 |
| 纯度 | Protein G purified |
| 克隆 | 单克隆 |
| 同种型 | IgG2a |
| 轻链类型 | kappa |

应用

Our [Abpromise guarantee](#) covers the use of **ab54673** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|-------|------|---|
| WB | | Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 34 kDa. |
| IHC-P | | Use a concentration of 5 µg/ml. |

靶标


| | |
|-------|---|
| 功能 | Likely to play a regulatory role in the development of the ventral forebrain. May play a role in craniofacial patterning and morphogenesis. |
| 序列相似性 | Belongs to the distal-less homeobox family. Contains 1 homeobox DNA-binding domain. |
| 细胞定位 | Nucleus. |

图片



Dlx2 antibody (ab54673) used in immunohistochemistry at 5ug/ml on formalin fixed and paraffin embedded human stomach.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Dlx2 antibody (ab54673)



Dlx2 antibody (ab54673) at 1ug/lane + NIH/3T3 cell lysate at 25ug/lane.

Western blot - Anti-Dlx2 antibody (ab54673)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors