

Product datasheet

Anti-Dlx2 antibody ab54673

2 图像

概述

产品名称	Anti-Dlx2抗体
描述	小鼠单克隆抗体to Dlx2
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Mouse, Human
免疫原	Recombinant fragment, corresponding to amino acids 1-110 of Human Dlx2

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab54673** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB

IHC-P

应用说明	WB: Use at a concentration of 1-5 µg/ml. IHC-P: Use at a concentration of 5 µg/ml.
------	---

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## 靶标

### 功能

Likely to play a regulatory role in the development of the ventral forebrain. May play a role in craniofacial patterning and morphogenesis.


### 序列相似性

Belongs to the distal-less homeobox family.  
Contains 1 homeobox DNA-binding domain.

### 细胞定位


Nucleus.

## 图片



IHC-P - Dlx2 antibody (ab54673)

Dlx2 antibody (ab54673) used in immunohistochemistry at 5ug/ml on formalin fixed and paraffin embedded human stomach.



Western blot - Dlx2 antibody (ab54673)

**Predicted band size : 34 kDa**

Dlx2 antibody (ab54673) at 1ug/lane +  
NIH/3T3 cell lysate at 25ug/lane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors