

Product datasheet

Anti-Filensin antibody ab54851

1 图像

概述

产品名称	Anti-Filensin抗体
描述	小鼠单克隆抗体to Filensin
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant fragment: EESRRPCAMV TPGAEEPSIP EPPKPAADQD GAEVLGTRSR SLPEKGPCKA LAYKTVEVVE SIEKISTESI QTYEETAVIV ETMIGKTKSD KKKSGEKS, corresponding to amino acids 567-665 of Human Filensin Run BLAST with ExPASy Run BLAST with NCBI

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab54851** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 75 kDa.

靶标

组织特异性	Lens.
疾病相关	Defects in BFSP1 are the cause of cataract cortical juvenile-onset (CCJO) [MIM:611391]. A juvenile-onset cataract with opacities restricted to the cortex of the lens, not involving the nucleus.
序列相似性	Belongs to the intermediate filament family.
细胞定位	Membrane. Cytoplasm. Cytoplasm > cytoskeleton. Membrane- and cytoskeleton-associated.

图片



Filensin antibody (ab54851) at 1 ug/lane + HL-60 cell lysate at 25ug/lane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors