

Product datasheet

Anti-MEIS2 antibody ab55647

2 References 2 图像

概述

产品名称	Anti-MEIS2抗体
描述	小鼠单克隆抗体to MEIS2
宿主	Mouse
经测试应用	适用于: WB, IHC-P, ICC/IF
种属反应性	与反应: Human
免疫原	Recombinant full length protein, corresponding to amino acids 1-382 of Human MEIS2

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG1
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab55647** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 52 kDa.
IHC-P		Use a concentration of 1 µg/ml.
ICC/IF		Use at an assay dependent concentration. PubMed: 21345011

## 靶标

### 组织特异性

Expressed in various tissues. In hematopoietic tissues, the lymphoid organs express high levels of MEIS2. Also expressed in some regions of the brain, such as the putamen.


### 序列相似性

Belongs to the TALE/MEIS homeobox family.  
Contains 1 homeobox DNA-binding domain.

### 细胞定位


Nucleus.

## 图片



IHC-P - MEIS2 antibody (ab55647)

MEIS2 antibody (ab55647) used in immunohistochemistry at 1ug/ml on formalin fixed and paraffin embedded human spleen tissue.



Western blot - MEIS2 antibody (ab55647)

MEIS2 antibody (ab55647) at 1ug/lane + HeLa cell lysate at 25ug/lane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors