

## Product datasheet

# Anti-SLC25A11 antibody ab56334

★★★★☆ 1 Abreviews 1 References 1 图像

### 概述

产品名称	Anti-SLC25A11抗体
描述	小鼠单克隆抗体to SLC25A11
宿主	Mouse
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Mouse, Recombinant fragment
免疫原	Recombinant full length protein, corresponding to amino acids 1-315 of Human SLC25A11

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG1
轻链类型	kappa

### 应用

Our [Abpromise guarantee](#) covers the use of **ab56334** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

应用	Ab评论	说明
----	------	----

IHC- P ★★★★★ 1/100. (see Abreview)


## 靶标

**功能** Catalyzes the transport of 2-oxoglutarate across the inner mitochondrial membrane in an electroneutral exchange for malate or other dicarboxylic acids, and plays an important role in several metabolic processes, including the malate-aspartate shuttle, the oxoglutarate/isocitrate shuttle, in gluconeogenesis from lactate, and in nitrogen metabolism.

**序列相似性** Belongs to the mitochondrial carrier family.  
Contains 3 Solcar repeats.

**细胞定位** Mitochondrion inner membrane.

## 图片



Western blot against tagged recombinant protein immunogen using ab56334 SLC25A11 antibody at 1ug/ml. Predicted band size of immunogen is 61 kDa

Western blot - SLC25A11 antibody (ab56334)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors