

Product datasheet

Anti-MLF2 antibody ab56378

2 图像

概述

产品名称	Anti-MLF2抗体
描述	小鼠单克隆抗体to MLF2
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant full length protein, corresponding to amino acids 1-249 of Human MLF2

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG1
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab56378** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB


应用说明 WB: Use at a concentration of 1-5 µg/ml.

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

靶标

组织特异性	Ubiquitously expressed.
序列相似性	Belongs to the MLF family.
细胞定位	Cytoplasm. Nucleus.

图片




Anti-MLF2 antibody (ab56378) at 5 µg/ml +
0.2 µg/lane of full length recombinant MLF2
protein with a tag

Secondary

Goat anti mouse IgG (H&L)-HRP conjugate at
1/5000 dilution

Predicted band size : 29 kDa

Observed band size : 53 kDa



All lanes : Anti-MLF2 antibody (ab56378) at
5 µg/ml

Lane 1 : 50 µg/lane of transfected 293T cell
line lysate

Lane 2 : 50 µg/lane of non-transfected 293T
cell line lysate

Secondary

Goat anti mouse IgG (H&L)-HRP conjugate at
1/2500 dilution

Predicted band size : 29 kDa

Observed band size : 30 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors