

Product datasheet

Anti-GRAP2 antibody ab56498

2 图像

概述

产品名称	Anti-GRAP2抗体
描述	小鼠单克隆抗体to GRAP2
宿主	Mouse
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human
免疫原	Recombinant fragment: HHFHQERRGG SLDINDGHCG TGLGSEMNAALMHRRTDPVQLQAAGRVRWARALYDFEAL EDDELGFHSG EVVEVLDSSN PSWWTGRLHN , corresponding to amino acids 226-316 of Human GRAP2 Run BLAST with ExPASy Run BLAST with NCBI

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab56498** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 38 kDa.

应用	Ab评论	说明
IHC-P		Use a concentration of 3 µg/ml.

靶标


功能	Interacts with SLP-76 to regulate NF-AT activation. Binds to tyrosine-phosphorylated shc.
序列相似性	Belongs to the GRB2/sem-5/DRK family. Contains 1 SH2 domain. Contains 2 SH3 domains.

图片



GRAP2 antibody (ab56498) used in immunohistochemistry at 3ug/ml on formalin fixed and paraffin embedded human spleen.

IHC-P - GRAP2 antibody (ab56498)



GRAP2 antibody (ab56498) at 1ug/lane + Jurkat cell lysate at 25ug/lane.

Western blot - GRAP2 antibody (ab56498)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors