

Product datasheet

Anti-GRAP2 antibody ab56498

2 图像

概述

产品名称	Anti-GRAP2抗体
描述	小鼠单克隆抗体to GRAP2
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human
免疫原	Recombinant fragment: HHFHQERRGG SLDINDGHCG TGLGSEMNAALMHRRHTDPVQLQAAGRVRW ARALYDFEAL EDDELGFHSG EVVEVLDSSN PSWWTGRLHN , corresponding to amino acids 226-316 of Human GRAP2 Run BLAST with ExPASy Run BLAST with NCBI

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab56498** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		
IHC-P		

应用说明

WB: Use at a concentration of 1-5 µg/ml.

IHC-P: Use at a concentration of 3 µg/ml.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

功能

Interacts with SLP-76 to regulate NF-AT activation. Binds to tyrosine-phosphorylated shc.


序列相似性

Belongs to the GRB2/sem-5/DRK family.

Contains 1 SH2 domain.


Contains 2 SH3 domains.

图片



IHC-P - GRAP2 antibody (ab56498)

GRAP2 antibody (ab56498) used in immunohistochemistry at 3µg/ml on formalin fixed and paraffin embedded human spleen.



Western blot - GRAP2 antibody (ab56498)

Predicted band size : 38 kDa

GRAP2 antibody (ab56498) at 1µg/lane + Jurkat cell lysate at 25µg/lane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors