

Product datasheet

Anti-KIF26A antibody ab57299

1 图像

概述

产品名称	Anti-KIF26A抗体
描述	小鼠单克隆抗体to KIF26A
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant full length protein, corresponding to amino acids 1-204 of Human KIF26A

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG1
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab57299** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB


应用说明 WB: Use at a concentration of 1-5 µg/ml.

Not yet tested in other applications.  
Optimal dilutions/concentrations should be determined by the end user.

## 靶标

<b>功能</b>	Atypical kinesin that plays a key role in enteric neuron development. Acts by repressing a cell growth signaling pathway in the enteric nervous system development, possibly via its interaction with GRB2 that prevents GRB2-binding to SHC, thereby attenuating the GDNF-Ret signaling. Binds to microtubules but lacks microtubule-based motility due to the absence of ATPase activity.
<b>序列相似性</b>	Belongs to the kinesin-like protein family. KIF26 subfamily. Contains 1 kinesin-motor domain.
<b>细胞定位</b>	Cytoplasm > cytoskeleton.

## 图片



**Predicted band size : 23 kDa**  
KIF26A antibody (ab57299) at 1ug/lane +  
HeLa cell lysate at 25ug/lane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors